

BioassayPlus Inc.

Windsor, ON, Phone: 313-455-0968, info@immunoassay.ca

Ready-Made Fly Food (Food Made Fresh Every Week)

Fly food B:

Similar to R, has more agar and cornmeal, less yeast, stores better.




Fly food M:

Similar to B, uses 20% more agar, food is firm.



Fly food R:

Uses Molasses as sugar source, soft, rich, suitable for weak stain.



Fly food A:

Similar to cornmeal, sucrose, yeast medium with hydroxyl benzoic acid methyl ester



Fly food J:

Similar to R, has more agar and cornmeal, less yeast, stores better.



YG food:



Banana Food:

Recipe by Cornell Stock Center



Grape Plates:

Recipe by Sarah Elgin. For feeding adults and collecting embryos.



Note: YG food recipe:

- 100g glucose/liter
- 15g Agar/liter 100g
- Brewer's Yeast/liter
- 3g Tygosept

Catalog #	Fly food type	\$ CAD
Food B		
8001-NV	Food B in narrow vials, 100vials/tray	\$49.00/tray
8001-WV	Food B in wide vials, 100vials/tray	\$59.00/tray
8001-BT	Food B in square bottom bottle with paper lid	\$2.20/bottle
8001-NVC	Food B in narrow vials, 100vials/tray with closure	\$55.00/tray
8001-WVC	Food B in wide vials, 100vials/tray with closure	\$66.00/tray
8001-NVYC	Food B in narrow vials, 100vials/tray with yeast and closure	\$59.00/tray
8001-WVYC	Food B in wide vials, 100vials/tray with yeast, and closure	\$69.00/tray
8001-BTYL	Food B in square bottom bottle with yeast, and paper lid	\$2.30/bottle

BioassayPlus Inc.

Windsor, ON, Phone: 313-455-0968, info@immunoassay.ca

8001-NV-PA	Food B with propionic acid only, in narrow vials, 100 vials/tray	\$49.00/tray
8001-WV-PA	Food B with propionic acid only in wide vials, 100vials/tray	\$59.00/tray
8001- BT-PA	Food B with propionic acid only in square bottom bottle	\$2.30/bottle
Food M		
8002-NV	Food M in narrow vials, 100vials/tray	\$53.00/tray
8002-WV	Food M in wide vials, 100vials/tray	\$64.00/tray
8002-BT	Food M in square bottom bottle	\$2.30/bottle
8002-NV-Teg	Food B with Tegosept only, in narrow vials, 100 vials/tray	\$53.00/tray
8002-BT-Teg	Food B with Tegosept only, in square bottom bottle	\$2.30/bottle
Food R		
8003-NV	Food R in narrow vials, 100vials/tray	\$49.00/tray
8003-WV	Food R in wide vials, 100vials/tray	\$59.00/tray
8003-BT	Food R in square bottom bottle	\$2.30/bottle
8003-NVC	Food R in narrow vials, 100vials/tray with closure	\$55.00/tray
8003-WVC	Food R in wide vials, 100vials/tray with closure	\$66.00/tray
8003-BT	Food R in square bottom bottle with paper lid	\$2.30/bottle
8004-NV-amp	Food R with ampicillin 100ug/ml, in narrow vials, 100 vials/tray	\$51.00/tray
8004-WV-amp	Food R with ampicillin 100ug/ml in wide vials, 100vials/tray	\$62.00/tray
8004-WV-amp	Food R with ampicillin 100ug/ml in wide vials, 100vials/tray with closure	\$68.00/tray
8004-BT-amp	Food R with ampicillin 100ug/ml in square bottom bottle	\$2.40/bottle
Food A		
8005-NV	Food A in narrow vials, 100 vials/tray	\$54.00/tray
8005-WV	Food A in wide vials, 100vials/tray	\$64.00/tray
8005-BT	Food A in square bottom bottle with paper lid	\$2.50/bottle
8005-PYNV	Food A in narrow vials, 100vials/tray with closure	\$59.00/tray
8005-PYWV	Food A in wide vials, 100vials/tray with closure	\$68.00/tray
Food J		
8006-NV	Food J in narrow vials, 100vials/tray	\$49.00/tray
8006-WV	Food J in wide vials, 100vials/tray	\$59.00/tray
8006-BT	Food J in square bottom bottle with paper lid	\$2.20/bottle
8006-PNV	Food J in narrow vials, 100vials/tray with closure	\$55.00/tray
8006-PWV	Food J in wide vials, 100vials/tray with closure	\$66.00/tray
8006-PBT	Food J in square bottom bottle with closure	\$2.20/bottle
Yeast Glucose food		
8007-NV	Yeast glucose food in narrow vials, 100vials/tray	\$49.00/tray
8007-WV	Yeast glucose food in wide vials, 100vials/tray	\$59.00/tray
8007-BT	Yeast glucose food in square bottom bottle	\$2.20/bottle
Banana Food		
8008-NV	Banana food in narrow vials, 100vials/tray	\$49.00/tray
8008-WV	Banana food food in wide vials, 100vials/tray	\$59.00/tray
8008-BT	Banana food food in square bottom bottle	\$2.20/bottle
Grape and Apple Juice Food		
8009-plate	Grape juice plate for feeding adult and collecting embyo 60 x 15mm	\$0.95/plate
8009-plate	Apple juice plate for feeding adult and collecting embyo 60 x 15mm	\$0.95/plate

Contact information: info@immunoassay.ca, Phone: 313-455-0968

Payment: please send PO to info@immunoassay.ca

Shipping: overnight shipment through Fedex

Free shipping for over 250.00 CAD order.